



Hi-Spec

FACTSHEET

Patient Lift 350kg Capacity



Patient lift showed in the fully stowed position.



Lift being deployed showing the toe flap and the trolley stop in the 'up' position and safety handrails.



'StepSafe' foot protection fitted as standard to allow the use of wanderlead controls.



Deployed with the toe ramp down and ready to receive a stretcher or carry chair.

Key Features

- Aluminium toe ramp and trolley stops anodised as standard and do not require painting.
- Mechanical override system in case of electrical failure.
- All key operations protected by electro-mechanical interlocks to prevent damage to the tail lift or the ambulance.
- The Patient Lift meets all UK and European CE Health & Safety requirements. StepSafe addresses BS EB 1756-2:2001 which is very specific with reference to 'crushing & shearing' particularly regarding danger to feet.

Britannia House | Junction Street | Darwen | Lancashire | BB3 2RB
Tel: 01254 700222 | Fax: 01254 700233 | Email: sales@hi-specnc.co.uk